FELINE LEUKEMIA VIRUS ANTIGEN-FELINE IMMUNODEFICIENCY VIRUS ANTIBODY TEST KIT (FeLV/FIV COMBO)

Feline Leukemia Virus Antigen-Feline Immunodeficiency Virus Antibody Test Kit (FeLV-FIV Combo) for the qualitative detection of feline leukemia virus Antigen and feline immunodeficiency virus antibodies in feline whole blood, serum, or plasma

INTENDED USE

The VETSCAN Feline Leukemia Virus Antigen-Feline Immunodeficiency Virus Antibody (FeLV-FIV Combo) Rapid Test is a visual rapid test for the qualitative detection of Feline Leukemia Virus antigen and Feline Immunodeficiency Virus antibody in feline whole blood, serum, or plasma. This test is for veterinary use only.

INSTRUCTIONS FOR USE

- Samples must be at room temperature (15-27 °C/59-80 °F), before running the assay – DO NOT HEAT.
- Previously frozen or older samples must be centrifuged before use.
- Serum or Plasma**, either fresh, previously frozen or stored at 2-7 °C (36-45 °F), may be used in this test. Serum or plasma may be stored for up to 7 days at 2-7 °C. For longer storage, sample should be frozen at -20 °C or colder (-4 °F).
- Whole Blood** may be used. Whole blood must be anticoagulated (e.g. EDTA, heparin) and may be used either fresh or after refrigeration at 2-7 °C (36-45 °F) for up to 5 days.
- Excessive hemolysis may obscure the results.
- EDTA or heparin in plasma will not affect the results.

PRECAUTIONS AND WARNINGS

- Important:** Do not remove Test Device from the pouch until ready for use.
- For veterinary use only.
- Do not use components after expiration date.
- The Test Device must be used as soon as possible after removing from pouch.
- The Test Device should be in a horizontal position on a flat surface while the test is performed.
- Use a separate Transfer Pipette for each test.
- The Chase Buffer is not interchangeable from serial (lot) to serial (lot).
- Contains sodium azide as a preservative.

STORAGE

- Test Devices and Chase Buffer must be stored at room temperature 15-27 °C (59-80 °F).
- Test Devices and Chase Buffer are stable until the expiration date when stored at 15-27 °C (59-80 °F).

KIT COMPONENTS

- Test Devices
- Chase Buffer Bottle
- Transfer Pipettes
- Instructions for Use

TEST PROCEDURE

- Remove the Test Device from the protective pouch and place on a flat surface. Label the Test Device with the subject I.D. or control identification. Only use the provided white space for this purpose. Do not cover up the test name.
- Gently mix the sample by inverting.
- Using a provided Transfer Pipette, hold the pipette vertically and dispense one drop of sample (whole blood, serum, or plasma) into sample well.
- Holding the Chase Buffer Bottle vertically, promptly add 2 drops of the chase buffer to the sample well.
- Read the results within 10 minutes. Do not read results after 10 minutes. Colored lines which appear after 10 minutes are not diagnostic and should be ignored.

INTERPRETATION OF TEST RESULTS

Positive results

A test is positive if a colored line appears under the FeLV (L) or FIV (I) areas and a line appears under the Control (C) area. Any intensity of the Test line L or I areas should be considered positive. Colored lines may be lighter or darker than each other.

Negative Results

The test is negative if only one line appears at the Control line (C) area.


Invalid Results


The test is invalid if no colored line appears at the Control line (C) area even if a colored line appears in the Test line L or I areas. If a colored line does not appear in the C area, the test is invalid and should be repeated.


ABAXIS CUSTOMER CARE


In Americas: +1 800 822 2947 . Outside Americas: +49 6155 780 210


Combo FeLV-FIV Test Procedure


- 

Add 1 drop of blood, serum or plasma to the sample well.
- 

Add 2 drops of Chase Buffer to the sample well.
- 

Read results within 10 minutes.
Positive Example
- 

Read results within 10 Minutes
Negative Example
- 

Read results within 10 minutes
Invalid Example
- 

Read results within 10 minutes.
Invalid Example

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Learn more about all of our products and services at www.abaxis.com Abaxis and VETSCAN are registered trademarks of Abaxis, Inc. ©Abaxis 2018.