Front

WITNESS[™] FeLV-FIV

FELINE LEUKEMIA VIRUS ANTIGEN – FELINE IMMUNODEFICIENCY VIRUS ANTIBODY TEST KIT

I. GENERAL INFORMATION
WITNESS[™] FeLV-FIV test is indicated for use when history and/or clinical signs may suggest an infectious feline retrovirus. It will also be particularly recommended prior to an FeLV vaccination, especially in cats belonging to an at-risk population.

II. TEST PRINCIPLES
The WITNESS[™] FeLV-FIV is a simple test based on Rapid Immuno Migration (RIM[®]) technology. FeLV antigen is detected using antibodies directed against the circulating p27 capsid protein. The detection of FIV is based on the presence of antibodies against the transmembrane region of the envelope protein along the strip. The complex is then captured on a sensitized nitrocellulose line where their accumulation causes the formation of a clearly visible pinkish band. Control bands located at the end of the reading window (C) ensure that the test was performed correctly.

III. SAMPLE COLLECTION
The test can be performed on unclotted whole blood anticoagulated with EDTA or heparin, serum or plasma.
• Hemolysis should be avoided – collect with a sterile needle and syringe.
• Hemolysis does not significantly interfere with the test, but strongly hemolyzed samples may partly obscure a weak positive band.

IV. SAMPLE STORAGE
Anticoagulated whole blood samples should preferably be tested immediately after collection but not longer than 4 hours after collection, if stored at room temperature.
• Reading is further delayed samples should be kept refrigerated at -20 °C and for up to 4 days.
• For Serum/Plasma: Samples (Serum and Plasma) must be kept frozen (20 °C).

V. KIT CONTENTS
• 10 to 25 Pouches, each containing 1 Test device and desiccant.
• 10 to 25 Pipettes.
• 1 to 2 Saline Buffer dropper bottles (5.4 mL).
• Instructions for use.

VI. PRECAUTIONS
Do not use this kit on any of its components after expiration date.
• Kit should be stored at +2 °C - 25 °C (32 °F - 77 °F). Kit should be frozen.
• Use the test immediately after opening the sealed pouch (within 10 minutes).
• Avoid touching or damaging the membrane in the sample wells or the results window.
• The WITNESS[™] device should be placed on a flat, horizontal surface while performing the test.
• Use a separate pipette for each sample.
• Hold pipette and buffer bottle vertically when dispensing sample and buffer respectively.
• Handle reagents and samples as biohazardous material.
• Saline buffer contains sodium azide as a preservative.
• For veterinary use only.

VII. TEST PROCEDURE
Important: Allow samples and buffer drops to fall onto the membrane in each sample window. Do not touch the pipette, sample or buffer drops, or buffer bottle top directly to the membrane.

I. SAMPLE APPLICATION
• Tear open a pouch provided and place the test device on a flat horizontal surface.
• Holding the pipette vertically, add one drop of sample to each test well (C).

II. BUFFER DISPENSING
• Check that the sample has fully absorbed into the membrane.
• Remove the cap from the buffer bottle, hold it vertically and add two drops of buffer to each sample well (C).
• Leave the test device for 10 minutes at room temperature. Sample/Buffer contact through the reading window.

3. READING TEST
After 10 minutes, observe the presence or absence of pinkish bands in the reading windows (C) and (D).
Note:
• It is possible to read the test before 10 minutes if two pinkish bands are clearly visible in the reading window (C) and (D).
• The presence of only one band in reading window (C) prior to the end of the development time (10 minutes), does not mean that the test is complete, as a test band may appear more slowly than the control band.

VIII. RESULTS
Valid Results
• Test is valid if a pinkish band is present in each reading window (C).

Interpretation of Results
• Positive: One band in reading window (C), with one band in window (D). Sample is positive for FeLV antigen and/or FIV antibodies.
• Negative: No band in reading window (C), with one band in window (D). Sample is negative for FeLV antigen and/or FIV antibodies.

Note:
• No band in control window (D) invalid test.
• A test result should always be interpreted in the context of all available clinical information and history for the cat being tested.

SYMBOL DESCRIPTIONS
Use by Date **ECREP** Authorized Representative in the European Community
Lot Code **LOT** Consult the instructions for use
SN Serial Number **SN** In vitro Diagnostics
Temperature Limitation **TL** Manufacturer

zoetis
Zoetis Inc.
Kalamazoo, MI 49007, USA
U.S. Vet. License No. 190
1-888-963-8401
www.zoetis.com

ECREP
SYMBIOTIS EUROPE SAS, a wholly-owned subsidiary of Zoetis
Genève Plaza, Bâle E.
23 Rue Pierre Gilles de Gennes, 69007
Lyon, FRANCE

WITNESS[™] FeLV-FIV

FELINE LEUKEMIA VIRUS ANTIGEN – FELINE IMMUNODEFICIENCY VIRUS ANTIBODY TEST KIT

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