

ZOETIS DIAGNOSTICS

vetscan VUE+

Automatically scan, display and share rapid test results with a next-generation reader

Eliminate manual processes and the margin for human error, and connect results directly to your PIMS via the Vetscan Hub™*, for a streamlined workflow that delivers reliable answers as well as peace of mind.



Reliable reader
**enables quick,
confident treatment
decisions**

- ✓ **Reliable results** you can trust
- ✓ **Reads Vetscan Rapid Tests**
- ✓ Part of the comprehensive **Zoetis Virtual Laboratory** ecosystem



Automated
process **increases
efficiency and frees
up staff time**

- ✓ Trusted turnkey solution — **insert a test and come back in 8-10 minutes†**
- ✓ **Timed and quick scan options** to suit your practice and patient needs
- ✓ Hardwired connection for **stable Hub, PIMS‡ and cloud integration**



Automation
and connectivity
with PIMS‡ **maximises
practice revenue**

- ✓ Transmits results to PIMS and **prevents missed charges**
- ✓ **Eliminates risk of misreads** and need for retesting
- ✓ Enables staff to **focus their valuable time elsewhere**

* Vetscan Hub required.

† Approximate. Refer to the instructions for the rapid test you are using for the specific read time.

‡ Compatible with select practice information management software (PIMS).


ZOETIS DIAGNOSTICS

vetscan Rapid™ Lateral Flow Tests

**Spend less time testing
and more time caring
for your patients**

With Vetscan® and
Witness® Rapid Tests

Get the flexibility and answers you need
to make fast, confident diagnoses with our
convenient lateral flow tests portfolio.



Reliable, point-of-care
rapid test **simplifies**
diagnoses

- ✓ Results in as little as **8 minutes**
- ✓ Delivered through **lateral flow technology**, proven for 30+ years



Easy operation
requires minimal
hands-on time

- ✓ **Simple, intuitive 3-step process** helps minimise errors
- ✓ **Automatic reading and reporting** with the Vetscan VUE+®†
- ✓ **Frees staff** to spend time where they are needed most








Efficient operation
and storage
minimises waste

- ✓ 12-month shelf life for most tests helps **ensure every test counts**
- ✓ Room temperature storage **saves valuable refrigerator space**

† Witness tests are not compatible with the Vetscan VUE+.

Choose from a wide portfolio of reliable rapid tests...

...and run only those needed to make the right clinical decisions.




Test	Detects antigens (Ag), antibodies (Ab) or hormones	Description	Time to results	Sensitivity	Specificity	Species
Canine Infectious Disease						
Witness Heartworm¹	Heartworm Ag	An efficient, easy-to-use visual rapid test for the detection of <i>Dirofilaria immitis</i> antigen in whole blood, serum or plasma.	10 minutes	97.0%	96.4%	
Vetscan Lepto²	<i>Leptospira</i> Ab	A simple, efficient rapid test to identify <i>leptospira</i> IgM antibodies as early as one week after infection.	10 minutes	83.7%	90.2%	
Feline Infectious Disease						
Witness FeLV/FIV^{3,4}	Feline Leukaemia Virus Ag, Feline Immunodeficiency Virus Ab	A simple, visual rapid test for the qualitative detection of FeLV antigens and FIV antibodies in whole blood, serum or plasma.	10 minutes	FeLV 92.9% FIV 93.8%	FeLV 96.5% FIV 93.4%	
Faecal						
Witness Parvo⁵	Parvovirus Ag	A simple in-clinic test for the qualitative detection of parvovirus antigens in canine feces.	10 minutes	86.3%	96.1%	
Witness Giardia⁶	Giardia intestinalis antigen in the faeces	A highly sensitive and specific test for the qualitative detection of giardia cyst antigens in canine faeces.	5 minutes	83%	100%	

Compatible with Vetscan VUE+ Lateral Test Flow Reader



Choose from a wide portfolio of reliable rapid tests...

...and run only those needed to make the right clinical decisions.

Test	Detects antigens (Ag), antibodies (Ab) or hormones	Description	Time to results	Sensitivity	Specificity	Species
Reproduction						
Vetscan Relaxin^{7,8}	Relaxin Hormone due to Pregnancy	An early, reliable point-of-care pregnancy test using serum or plasma.	10 minutes	Canine: 99% Feline: 100%	Canine: 100% Feline: 95.9%	
Vetscan LH⁹	Canine and Feline Luteinising Hormone	A simple in-clinic serum rapid test to identify the luteinising hormone (LH) for optimal breeding that may also help distinguish between ovariectomised females from those that are sexually intact.	20 minutes	Feline: 89.3%	Feline: 92.6%	
Test	Detects	Description	Time to results	Level of Agreement	Comparison Standard	Species
Pancreatic						
Vetscan cPL^{5,10}	Canine Pancreatic Lipase	A simple and efficient in-clinic test for the semi-quantitative detection of pancreas-specific lipase in canine serum or plasma.	10 minutes	ICC*: 0.96	Spec cPL	

⁹ Must be read with the Vetscan Vue Rapid Test Analyser

¹⁰ May also aid in distinguishing between ovariectomised female dogs or cats from those that are sexually intact

* Intraclass Correlation Coefficient (Agreement Scoring)

Compatible with Vetcan VUE+ Lateral Test Flow Reader



Data from USDA submission or recent peer reviewed publications.

References: 1. Starkey, LA, et al., Parasit Vectors, 2017 2. Study Report No. D866R-US-14-015-Supplement (D6015Q005) 3. Outcomes Research Report No. D880Z-US-14-014 4. Outcomes Research Report No. D880Z-US-14-013 5. Decaro, N, et al., Vet J, 2013. 6. Zoetis data on file: Giardia Antigen Test study, Tierärztliches Laboratorium, Freiburg, Germany, February 2013. Cyst observation under microscope was used as reference method. 100 positive dog samples and 50 negative dog samples were used as reference sample panel. 7. DiGangi, BA, et al., J Am Vet Med Assoc, 2010 8. Study Report: Witness Relaxin Global Dossier Version n°1 9. Krecic MR, DiGangi BA, Griffin B. Accuracy of a point-of-care luteinizing hormone test for help in distinguishing between sexually intact and ovariectomized or castrated domestic cats. Journal of Feline Medicine and Surgery. 2018;20(10):955-961. 10. Cridge, H, et al., J Vet Intern Med, 2018

Zoetis Australia Pty Ltd. ABN 94 156 476 425. Level 6, 5 Rider Boulevard, Rhodes, NSW 2138. © 2025 Zoetis Inc. All rights reserved. April 2025. MM-40696

