

Routine Cleaning Guide




vetscan IMAGYST™

LOOK DEEPER

zoetis

Routine Cleaning

1. **Turn off** the scanner by pressing the button
2. Remove the **power cable** and the **network cable** from the scanner
3. **Slide** the **topmost plate** to the **front**




4. **Wipe** all **surfaces** with a **soft, damp, gauze swab and distilled water**.
Use cotton swab with distilled water for cleaning the edges
5. **Dry** all surfaces with a **soft, damp, gauze swab or Kimwipes™**
6. **Slide** the **topmost plate** back to its normal position
7. Place a **protective cloth** on the **glass** beneath the **objective**
8. You can first try to clean the objective in place, using a microfiber cloth.
If necessary, try adding **warm distilled water** to the **cloth or using a cotton swab**

IMPORTANT: Never pour or spray any liquids directly on the scanner

9. **Connect** the **cables** and **switch the scanner on**
10. Verify by completing the scan again and validating the images
If the scans are still blurry follow the **In-Depth Cleaning** procedure

In-Depth Cleaning

1. **Turn off** the scanner by pressing the button
2. Remove the **power cable** and the **network cable** from the scanner
3. **Slide** the **topmost plate** to the **front**
4. **Remove** the **black casing on the cover plate** and wipe clean with a **soft damp, gauze swab**
5. **Wipe** all surfaces with a **soft, damp, gauze swab**. You can use a **microscope cleaner**. Use cotton swab with distilled water for cleaning the edges
6. **Replace** the **top plate cover**
7. **Slide** the **topmost plate** back to its normal position
8. **Connect** the cables and **switch the scanner on**
9. **Verify by completing the scan again** and validating the images:
 - a. *The initial scan after cleaning produces a good clear image*
 - b. *The live view from the camera produces good quality image*
 - c. *Scanning of a know sample slide produces good quality image*
10. If the quality is still poor, redo all steps and use **Isopropyl Alcohol** for the cleaning and follow the In-Depth Cleaning procedure for the objective




In-Depth Cleaning **the objective**

Should the objective become dirty as a result of faecal sample material or immersion oil, perform these steps to clean it:

1. **Turn off** the scanner by pressing the button
2. Remove the **power cable** and the **network cable** from the scanner
3. Place a **protective cloth** on the **glass** beneath the **objective***
4. You can first try to clean the objective in place, using a microfiber cloth.
If necessary, try adding **warm distilled water** to the **cloth** or using a **cotton swab**

IMPORTANT: Never pour or spray any liquids directly on the scanner

5. If the **objective** does not come clean, carefully unscrew it



6. **Use the microfibre cloth and warm distilled water to clean the lens**

1. **IMPORTANT:** To clean immersion oil from the lens, use warm distilled water, isopropyl alcohol or suitable microscope cleaning fluid. Use isopropyl alcohol or wipes sparingly
2. **IMPORTANT:** Using **acetone or xylene** to clean the lens may **damage it**; however, if the lens is covered with glue/adhesive, cleaning the lens may require using stronger cleaners*

7. Use a **bulb syringe** to remove any dust from the lens and scanner
8. Replace the objective by gently screwing it back in place

IMPORTANT: Do not forcefully screw the objective back onto the scanner

9. Remove the cloth covering the glass, replace the power cable, replace the network cable, and **turn the scanner on**
10. Check the scanner by running a **test scan**

*For further help or If you are unsure about what cleaner to use, please reach out to your Zoetis Diagnostic Technical Support team to clarify

Further information can be obtained from Zoetis UK Ltd. Birchwood Building, Springfield Drive, Leatherhead, Surrey, KT22 7LP.
www.zoetis.co.uk • Customer support 0345 300 8034 • CustomerSupportUK@zoetis.com • Produced April 2024 • MM 31383

zoetis