

vetscan OptiCell

Accurate automated complete blood count results with advanced AI technology¹

A validation study evaluated the accuracy and precision* of Vetscan OptiCell™ to the Advia® 2120 Reference Laboratory Hematology Analyzer, according to ASVCP guidelines.



Accuracy Study¹



The study analyzed anticoagulated (EDTA) blood samples



COLLECTED FROM
113 canines
72 felines



TESTED WITH
2 Vetscan OptiCell
analyzers with
randomized operators



TESTED FOR
Method comparison with the
Advia 2120 and 200-cell
manual blood count by a clinical
pathologist (for subset of samples)

Study Results¹

Vetscan Opticell demonstrated excellent agreement for key parameters[†] vs. the Advia 2120.

Accuracy measurements were determined using Lin's concordance correlation coefficient for most parameters (ρ):^{2,3}

INTERPRETATION OF LEVEL OF AGREEMENT³

- Excellent agreement beyond chance**
>0.75
- Good agreement beyond chance**
0.40-0.75
- Poor agreement beyond chance**
<0.40

MINIMAL BIAS




Canine
Minor proportional and constant
bias for key analytes




Feline
Minor proportional and constant
bias for key analytes


WBC ($10^9/L$)




RETIC ($10^3/L$)



RBC ($10^{12}/L$)



PLT ($10^9/L$)



WBC=White Blood Cells, RETIC=Reticulocytes, RBC=Red Blood Cells, PLT=Platelets, ASVCP=American Society of Veterinary Clinical Pathology.

* High accuracy is achieved when an instrument hits the correct target or value. High precision is achieved when an instrument generates consistent results from repeated tests.

[†] Key parameters include: WBC, RBC, PLT and RETIC.

Reference: 1. Data on file, Study No. DHXZ-US-24-235, 2024, Zoetis, Inc. 2. Dawson B, Trapp RG. Basic and clinical biostatistics. 4th ed. Lange Medical Books; 2004:159-174. 3. Landis JR, Koch GG. The measurement of observer agreement for categorical data. Biometrics. 1977;33(1). <https://doi.org/10.2307/2529310>
<https://doi.org/10.2307/2529310>

This document contains examples of performance data and is not exhaustive of analysis that was performed in the validation study. Additional data available upon request.

Powerful AI and Viscoelastic Focusing technology enable precise automated complete blood count results¹

Precision Study



The study analyzed canine EDTA blood samples and tri-level commercial quality control material



TESTED FOR
Precision with canine samples
Precision with tri-level control material



TESTED WITH
2 Vetscan OptiCell analyzers and 40 replicates across both analyzers using 3 control levels



TESTED ON
148 repeat runs using canine samples

Study Results

Vetscan OptiCell demonstrated precise results for key parameters, including PLT, across control levels and clinical samples.¹

<8.5%

Total precision below 8.5% for most parameters

<5%

Replicate precision was below 5% for most parameters

<6.5%

Clinical precision below 6.5% for most parameters

Short term precision with controls¹

Parameter	Control Level	Mean	SD Total	CV% (Replicate)
RBC (M/ μ L)	Low	2.67	0.06	2.3
WBC (K/ μ L)	High	19.41	0.97	3.8
PLT (K/ μ L)	Low	70.67	2.87	3.5
NEU (K/ μ L)	High	10.43	0.53	4.5
LYM (K/ μ L)	High	7.04	0.40	3.8
MON (K/ μ L)	High	1.43	0.14	7.6

Wide reportable ranges across all parameters¹

HGB (g/dl): Analyzer A+B



Reportable range analyzed for RBC, WBC, HCT, HGB, and RETICs. Only HGB shown here for example.

Vetscan OptiCell delivers advanced cell focusing with AI-powered technology that evaluates blood cells in their native state

Reference-lab quality

Automated complete blood count results¹

Viscoelastic Focusing

Enables a cartridge-based design to minimize errors

Detailed flags

Identify abnormal cell morphology

HGB=Hemoglobin, NEU=Neutrophils, LYM=Lymphocytes, MON=Monocytes, SD=standard deviation, CV=coefficient of variation.

Reference: 1. Data on file, Study No. DHXNZ-US-24-235, 2024, Zoetis, Inc.

This document contains examples of performance data and is not exhaustive of analysis that was performed in the validation study. Additional data available upon request.

LOOK DEEPER

Learn more at zoetisdiagnostics.com/opticell

All trademarks are the property of Zoetis Services LLC or a related company or a licensor unless otherwise noted. Advia is a registered trademark of Siemens Healthcare Diagnostics Inc. © 2024 Zoetis Services LLC. All rights reserved. VTS-01553

zoetis