

Diabetes in pets: clinical signs and diagnosis

Diabetes mellitus occurs when insulin is absent or not functioning properly, leading to a pet not being able to use glucose normally. Signs include:

- Excessive drinking
- Excessive urination
- Increased appetite
- Weight loss
- Lack of energy or depression
- Weakness of the back legs

Lifestyle changes to help your pet

Together with your veterinarian you can help manage your pet's diabetes with proper schedules and health tracking.

- Work with your vet to establish a dietary and exercise regime to help regulate blood sugar levels
- Track your pet's glucose and weight regularly, informing your vet of any changes

Accurate blood glucose results trusted by veterinarians^{1,2}

Find peace of mind knowing you're providing the best care for your cat or dog with an at-home testing monitor that's trusted and used by vets.



AlphaTrak 3 includes:

- Reliable test results your vet can depend on for management
- A bluetooth-connected app that makes it easy to track, review and graph your pet's glucose levels and share them with your vet
- Support, including instruction videos and one-on-one live help from Zoetis



References

1. TI-08743, Study report: D876R-US-21-040
2. AlphaTrak 3 Claim Statements Research, MedSurvey and Zoetis US Market Research 2022.

Zoetis Australia Pty Ltd. ABN 94 156 476 425. Level 6, 5 Rider Boulevard, Rhodes NSW 2138. © 2021 Zoetis Inc. All rights reserved. December 2022 MM-23794

zoetis AlphaTrak3

AlphaTrak[®]3 Trusted support for your diabetic pet

Partner with your veterinarian to confidently manage your cat or dog's diabetes through regular glucose monitoring, healthy habits, and prescribed treatments.



Improve care with at-home monitoring

Monitoring your pet's blood glucose helps track abnormal changes, and gives your veterinarian necessary information for proper treatment.

Benefits of at-home glucose monitoring

- Supports your vet with information to manage your pet's treatment plan
- Helps lower the likelihood of in-clinic glucose fluctuations due to stress
- Can reduce additional trips to the vet

Proper dosing for effective treatment

- Closely follow your vet's instructions on the specific insulin product you should use, syringe type, and dose
- Keep a log of the amount of insulin and time when each injection is given



Get an accurate picture with AlphaTrak³

Pet diabetes is similar to human diabetes, with a key difference in how glucose is distributed in the blood. Therefore, it is vital to use a home glucose monitor that is calibrated for pets to get the most accurate readings.

AlphaTrak 3 is designed specifically for dogs and cats to ensure accurate test results.¹

Simple testing for you and your pet



Auto on

Test strips turn the meter on and off for convenient testing



Small sample size

Only one small drop of blood needed (0.3 uL)



Quick testing

No wait time, results in seconds



Saved automatically

Easily track data in the meter



One-handed disposal

Patented eject button for simple lancet removal



Track insulin management

A digital log that's easy to access



Easily save and share data

Use the AlphaTrak 3 app to automatically track, review and graph glucose results, record your pet's daily activities, set up reminders and share data with your vet.

