

# AlphaTrak<sup>®</sup>3

## Accurate blood glucose monitoring for optimal care<sup>1,2</sup>

AlphaTrak 3 makes testing and monitoring blood glucose easier for everyone.

**AlphaTrak 3 is specifically calibrated for dogs, cats, and horses**

Humans, dogs, cats, and horses each have differences in the distribution of glucose in plasma and red blood cells.

Using a human meter in a dog, cat or horse can lead to an underestimation of blood glucose concentrations and inaccurate readings.<sup>3</sup>



 **Bluetooth<sup>®</sup>**



**Access accurate results with an extended blood glucose testing range<sup>1</sup>**

AlphaTrak 3 has a testing range of 1.1-41.7 mmol/L to allow for more pinpointed results for cats, dogs, and horses.



**zoetis**



# Testing made simple

Easily monitor blood glucose levels of dogs, cats and horses in the clinic, stable-side or from the comfort of your own home.



## Auto on

Test strips automatically turn the meter on and off for convenient testing



## Small sample size

Only 0.3 uL of blood needed



## Quick testing

No wait time, results in seconds



## Saved automatically to app

Digital data reduces the risk of human error or missed results, and makes it easy to trend and share

# Efficient patient tracking

Save time at your practice with enhanced digital results and tracking support. Connecting with client data allows for improved treatment and greater client peace of mind.



## Seamless data capture

Save glucose results and build curves easily




## Store results in PIMS

Connects to most practice management software via Wi-Fi with FUSE/HUB




## Easy data transfer

Request clients to share glucose data or reports from at-home testing as needed to incorporate into the patient's medical record



Desktop web app for your practice



Mobile app for clients and practice



# Improve patient outcomes

See how AlphaTrak 3 can help you partner with clients to increase the quality of care for your patients

Learn more at [zoetisdiagnostics.com/ie](https://zoetisdiagnostics.com/ie)  
Contact Zoetis: [VetScanIreland@Zoetis.com](mailto:VetScanIreland@Zoetis.com)

## References:

1. Data on file, Study Report No. D876RUS-21-040, Zoetis 2022.
2. Data on file, Study Report D850R-US-22-005, Zoetis Inc.
3. Data on file, Study Report No. 11-06-MC-D-CT-AT, Zoetis 2022.

