

# AlphaTrak2 vs AlphaTrak3

## Features & Benefits



**AlphaTRAK**

AlphaTrak3



Within 60 seconds	Blood Collection Time	No limits	Less strips required / wasted
Double Sided	Testing Strip Design	One Sided	Easier to use - more obvious where to put the sample
0.3 µL	Sample Size	0.3 µL	
10 seconds	Results Time	5 seconds	Much quicker, especially if you are under a time pressure
20-750 mg per dL	Glucose Measurement Range	20-750 mg per dL	
	App & Bluetooth Connectivity	Yes	View your patients results easily:
	Data Sharing Capabilities	Yes	<ul style="list-style-type: none"> <li>- New pet owner &amp; Vet app</li> <li>- Results transferred via Bluetooth</li> <li>- Automatic blood glucose curve created</li> </ul>



Introducing the **new** AlphaTrak3 blood glucose monitoring system from Zoetis



Further information can be obtained from Zoetis UK Ltd. Birchwood Building, Springfield Drive, Leatherhead, Surrey, KT22 7LP.  
www.zoetis.co.uk • Customer support 0345 300 8034 • CustomerSupportUK@zoetis.com • Produced April 2023 • MM-25404

**zoetis**