




REMEMBER THEY NEED TAKING CARE OF TOO

Scan now to go to the
VETSCAN YouTube channel
for help with your regular
maintenance checks.
youtube.com/vetscan



HM5 Maintenance Record

Serial Number:

Year:

Weekly Maintenance | Wash Head Cleaning

1	2	3	4	5	6	7	8	9	10	11	12	13
14	15	16	17	18	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	36	38	39
40	41	42	43	44	45	46	47	48	49	50	51	52

Bi-weekly Maintenance | Soak Cleaning with HemaClean*

January		February		March		April		May		June		July		August		September		October		November		December	
Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3	Week 1	Week 3

It is recommended a Quality Control is run every time you change a reagent pack

VS2 Maintenance Record

Serial Number:

Year:

Weekly Maintenance | Foam Filter Cleaning

1	2	3	4	5	6	7	8	9	10	11	12	13
14	15	16	17	18	19	20	21	22	23	24	25	26
27	28	29	30	31	32	33	34	35	36	36	38	39
40	41	42	43	44	45	46	47	48	49	50	51	52

For Technical Support please contact 0345 300 8034

*Wear gloves when handling HemaClean solution and Quality Control vial.

Further information can be obtained from Zoetis UK Limited, 1st Floor, Birchwood Building, Springfield Drive, Leatherhead, Surrey, KT22 7LP. www.zoetis.co.uk Customer support: 0345 300 8034 or CustomerSupportUK@zoetis.com
Date of preparation: June 2022. MM-12764